

Request for IH sampling media equipment

Industrial Hygiene Laboratory Services from Liberty Mutual Insurance

Ship attention

E-Mail

Phone number

Order date

Date needed*

Company

Street address

City

State

Zip

Sampling media

Please indicate quantities needed and include **blanks** in your total quantity.

Qty Sorbent Tubes

_____	Charcoal, small (SKC 226-01)
_____	Charcoal, large (SKC 226-09)
_____	ORBO 34 (hydrogen sulfide)
_____	<input type="text"/> (ammonia – choose media)
_____	ORBO 78 (ETO)
_____	Sep paks, WAT047205 (aldehydes) ^{1,2}
_____	SKC 226-10 (amines)
_____	SKC 226-10-03 (alkanolamines)
_____	SKC 226-10-03 & 25mm MCE prefilter (inorganic acids)
_____	SKC 226-10-03 & PTFE prefilter (organic acids)
_____	SKC 226-119, DNPH (aldehydes) ^{1,2}
_____	SKC 226-17-1A, small carulite (mercury)
_____	SKC 226-17-3A, large carulite (mercury)
_____	SKC 226-183, order by set: 2 tubes in series (diacetyl) ⁴
_____	SKC 226-193-UC, use w/ 225-9030 (peracetic acid)
_____	SKC 226-210, Soda-lime (hydrogen cyanide)
_____	SKC 226-30-04 with PTFE Filter & foil (PNAs)
_____	SKC 226-30-08, order by set: 2 tubes in series (acrylic acid)
_____	SKC 226-30-18, NITC (EDA, DETA, TETA)
_____	SKC 226-39 (nicotine)
_____	SKC 226-39 with Swinnex Filter (PCB)
_____	SKC 226-40-02 (nitrogen dioxide)
_____	SKC 226-40A, order by set: 3 tubes in series (NO/NO ₂) ¹
_____	SKC 226-57 (ethylene glycol) ²
_____	SKC 226-80, IABC (sulfur dioxide) ²
_____	SKC 226-82, order by set: 2 tubes in series (methanol)
_____	SKC 226-95, XAD 7 (phenol)
_____	SKC 226-98, acid treated XAD 7 (cyanoacrylates)

Other:

Notes

- Media must be kept cold
 - Limited shelf-life
 - Please indicate if you need the filters pre-weighed
 - Media must be wrapped in aluminum foil after sampling
 - Media charges apply
- * UPS Ground shipping is free.

Wipe sampling kits

 Select from pull down menu:

(Indicate number of wipes, templates, gloves needed)

Qty Common Filters

_____	MCE: 25mm (asbestos)		
_____	MCE: 2-pc, 37mm (metals)		
_____	PVC 25mm, 2-pc pre-weighed		
_____	PVC 2-pc, 37mm pre-weighed		
_____	PVC 3-pc, 37mm pre-weighed		
_____	PPI w/PVC ^{3,5}	Preweighted?	<input type="text"/>
_____	IOM-PVC ^{3,5}	Preweighted?	<input type="text"/>
_____		Calibrator?	<input type="text"/>

Qty Specialty Filters

_____	GFF with foil (coal tar pitch volatiles)		
_____	GFF, NaNO ₂ with foil (ozone) ^{2,4}		
_____	GFF with 1,2 p,p (isocyanates) ^{1,2}		
_____	IOM-GFF (resin acids) ⁵	Calibrator?	<input type="text"/>
_____	IOM-PTFE (pharma) ⁵	Calibrator?	<input type="text"/>
_____	PPI-MCE Thoracic (sulfuric acid) ⁵		
_____	PTFE, 37mm (metal working fluid/asphalt fume)		
_____	PTFE/Silver membrane (chlorine & bromine)		
_____	SKC 225-9030 w/tin foil (hydrogen peroxide) ^{4,5}		
_____	SKC 225-9031 (fluorides, aerosol & gases) ^{1,2,5}		
_____	SKC 225-9032 (hydrogen chloride/bromide, nitric) ^{2,5}		
_____	SKC 225-9033 (phosphoric, sulfuric) ^{2,5}		
_____	SKC 225-317 (diesel particulate matter)		

Qty Passive Vapor Monitors

_____	3M 3500 badge ²
_____	3M 3520 badge ²
_____	3M 3551 badge (ETO) ²
_____	N555AT, Assay Technology badge (ETO) ²
_____	N566AT, Assay Technology badge (solvents) ²
_____	N571AT, Assay badges (aldehydes) ¹
_____	X575AT, Assay badges (nitrous oxide) ¹

Miscellaneous

Qty

_____	Cooler w/ ice packs
_____	Fritted Bubblers
_____	KI solution for chlorine dioxide

Qty

_____	Glass Vials w/ parafilm (liquid bulks)
_____	Plastic Vials (solid bulks)
_____	Vials for isocyanate field extraction

Equipment

Requires completed Equipment terms and conditions, and request form (RC6128) prior to shipment

Gas testing

Choose type: Breathing air Pipeline

- _____ Cylinders
- _____ 47mm teflon filters
- _____ Fixtures

Pump loan

(Hoses, clips, and chargers included)

- _____ Pumps rate (Lpm): _____
- _____ Flow-limiters
- _____ cyclones (use 2pc cassette, rate: 1.7 Lpm)

Comments/other: (Qty, Description, Sampling for...)

Liberty Mutual Industrial Hygiene Lab
94 South Street, Hopkinton MA 01748
Telephone: (800)230-6263 Fax: (508) 435-3575
E-Mail this form to: LMIHLaboratory@LibertyMutual.com